

## Product Datasheet

### HSD17B4 peptide (orb234507)

<b>Description</b>	The peptide is used to anti-HSD17B4 antibody reactivity
<b>Reactivity</b>	Human, Mouse, Primate, Rat
<b>Tested Applications</b>	Blocking
<b>Target</b>	HSD17B4
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb234876 (HSD17B4 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">P51660</a> , <a href="#">P51659</a> , <a href="#">P97852</a>
<b>Entrez</b>	<a href="#">15488</a> , <a href="#">79244</a> , <a href="#">3295</a>
<b>Expiration Date</b>	12 months from date of receipt.