

Product Datasheet

NUP210L peptide (orb234559)

Description	The peptide is used to anti-NUP210L antibody reactivity
Reactivity	Human, Rat
Tested Applications	Blocking
Target	NUP210L
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb234928 (NUP210L antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q5VU65 , Q9D2F7
Entrez	91181 , 77595
Expiration Date	12 months from date of receipt.