

Product Datasheet

ACC alpha (pS80) peptide (orb234621)

Description	The peptide is used to anti-ACC alpha (pS80) antibody reactivity
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	ACACA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb234990 (ACC alpha (phospho-S80) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q5SWU9 , Q13085 , P11497
Entrez	60581 , 107476 , 31
Expiration Date	12 months from date of receipt.