

Product Datasheet

CaMK2 alpha/delta (pT286) peptide (orb234687)

Description	The peptide is used to anti-CaMK2 alpha/delta (pT286) antibody reactivity
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	CAMK2A; CAMK2D
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb235056 (CaMK2 alpha/delta (phospho-T286) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q13557 , P15791 , P11798 , Q9UQM7 , P11275 , Q6PHZ2
Entrez	108058 , 12322 , 25400 , 815 , 24246 , 817
Expiration Date	12 months from date of receipt.