

## Product Datasheet

### RPA2 Antibody (PE/Cy7) (orb2361683)

<b>Description</b>	RPA2 Antibody (PE/Cy7)
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	PE/Cy7
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Recombinant human RPA2 Protein
<b>Target</b>	RPA2
<b>Preservatives</b>	PBS with 0.02% Proclin300.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1,K
<b>Clonality</b>	Monoclonal
<b>MW</b>	29 kDa
<b>Uniprot ID</b>	<a href="#">P15927</a>
<b>Dilution Range</b>	WB=1:500-2000
<b>Expiration Date</b>	12 months from date of receipt.