

Product Datasheet

RPA2 Antibody (PE/Cy5.5) (orb2361684)

Description	RPA2 Antibody (PE/Cy5.5)
Species/Host	Mouse
Conjugation	PE/Cy5.5
Tested Applications	WB
Immunogen	Recombinant human RPA2 Protein
Target	RPA2
Preservatives	PBS with 0.02% Proclin300.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1,K
Clonality	Monoclonal
MW	29 kDa
Uniprot ID	P15927
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.