



Product Datasheet

AKT1 Antibody (orb238167)


Description	AKT1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human RAC-alpha serine/threonine-protein kinase protein (1-480AA)
Target	AKT1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31749
Dilution Range	WB:1:500-1:5000, IHC-P:1:100-1:1000, IF:1:200-1:500, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Colo320(Lane 1)...



Immunohistochemical staining of human pr...



Western blot analysis of human tonsil ti...