

Product Datasheet

Death domain-containing protein CRADD antibody (orb238226)

Description	Death domain-containing protein CRADD antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Death domain-containing protein CRADD protein (1-199AA)
Target	CRADD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P78560
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa(Lane 1), K...