
Product Datasheet


**Dolichyl-diphosphooligosaccharide-protein
glycosyltransferase subunit DAD1 antibody
(orb238234)**

Description

Dolichyl-diphosphooligosaccharide-protein



antibody


Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated


Tested Applications

ELISA, IF, IHC

Western blot analysis of EC109(Lane 1), ...

Immunogen

Recombinant Human Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit DAD1 protein (1-113AA)


Target

DAD1

Preservatives

 Preservative: 0.03% Proclin 300.
 Constituents: 50% Glycerol

Immunohistochemical staining of human pr...

Form/Appearance

Liquid

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality


Polyclonal

Uniprot ID
P61803
Dilution Range

IHC-P:1:20-1:200, IF:1:50-1:200

Expiration Date

12 months from date of receipt.



Immunofluorescence analysis of HepG2 cel...