

Product Datasheet

Dolichyl-diphosphooligosaccharide-protein glycosyltransferase subunit DAD1 antibody (Biotin) (orb238235)

Description	Dolichyl-diphosphooligosaccharide-protein glycosyltransferase subunit DAD1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit DAD1 protein (1-113AA)
Target	DAD1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P61803
Expiration Date	12 months from date of receipt.