

Product Datasheet

Glyceraldehyde-3-phosphate dehydrogenase antibody (HRP) (orb238249)

Description	Glyceraldehyde-3-phosphate dehydrogenase antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Glyceraldehyde-3-phosphate dehydrogenase protein (3-335AA)
Target	Ovomucoid
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04406
Dilution Range	WB:1:1000-5000, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.