

Product Datasheet

Programmed cell death protein 6 antibody (orb238258)

Description Programmed cell death protein 6 antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC

Immunogen Recombinant Human Programmed cell death protein 6 protein (1-191AA)

Target PDCD6

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Note For research use only


Isotype IgG

Clonality Polyclonal


Uniprot ID [O75340](#)

Dilution Range IHC-P:1:200-1:500, IF:1:50-1:200


Expiration Date 12 months from date of receipt.



Western blot analysis of EC109(Lane 1), ...



Immunofluorescence analysis of PC-3 cell...



Immunohistochemical staining of human lu...