



Product Datasheet

threonine-protein kinase PAK 4 antibody (orb238269)


Description	threonine-protein kinase PAK 4 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Serine/threonine-protein kinase PAK 4 protein (1-426AA)
Target	PAK4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O96013
Dilution Range	WB:1:500-1:5000, IHC-P:1:500-1:1000, IF:1:200-1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human prostate tissue.



Immunofluorescence analysis of Hela cells.



Immunohistochemical staining of human brain tissue.