



Product Datasheet

Tyrosine-protein phosphatase non-R type 6 antibody (orb238461)


Description	Tyrosine-protein phosphatase non-R type 6 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Tyrosine-protein phosphatase non-receptor type 6 protein (1-597AA of Isoform 2)
Target	PTPN6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P29350
Dilution Range	WB:1:1000-1:5000, IHC-P:1:20-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat(Lane 1),...



Immunohistochemical staining of human to...



Immunohistochemical staining of human sp...