

Product Datasheet

Phospho-MYL12B (Thr19 + Ser20) Rabbit Polyclonal Antibody (APC-Cy5.5) (orb2384899)

Description	Phospho-MYL12B (Thr19 + Ser20) Rabbit Polyclonal Antibody (APC-Cy5.5)
Species/Host	Rabbit
Conjugation	APC/Cy5.5
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Mylc2b around the phosphorylation site of Thr19 + Ser20 RA(p-T)(p-S)N
Target	MYL12B
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
Dilution Range	ICC/IF=1:50-200, IF=1:100-500
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)