



## Product Datasheet

### alpha-Galactosidase A antibody (Biotin) (orb238502)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | alpha-Galactosidase A antibody (Biotin)                          |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human  |
| <b>Conjugation</b>         | Biotin   |
| <b>Tested Applications</b> | ELISA  |
| <b>Immunogen</b>           | Recombinant Human Alpha-galactosidase A protein (32-429AA)       |
| <b>Target</b>              | GLA  |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol   |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C.<br>Avoid repeated freeze. |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <b>P06280</b>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                  |



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...