



Product Datasheet

Annexin A5 antibody (orb238513)


| | |
|----------------------------|---|
| Description | Annexin A5 antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC, WB |
| Immunogen | Recombinant Human Annexin A5 protein (2-320AA) |
| Target | ANXA5 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P08758 |
| Dilution Range | WB:1:1000-1:5000, IHC-P:1:20-1:200 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of EC109(Lane 1), ...



Western blot analysis of HepG-2 whole ce...



Immunohistochemical staining of human st...