



Product Datasheet

Carbonic anhydrase 1 antibody (orb238521)


Description	Carbonic anhydrase 1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human Carbonic anhydrase 1 protein (2-261AA)
Target	CA1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P00915
Dilution Range	WB:1:500-1:5000, IHC-P:1:20-1:200, IF:1:50-1:200, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1):K562 w...



Immunoprecipitation analysis of K562 who...



Western blot analysis of K562 whole cell...