



Product Datasheet

Stromelysin-2 antibody (orb238978)


Description	Stromelysin-2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Stromelysin-2 protein (99-476AA)
Target	MMP10
Preservatives	Preservative: 0.1mol/L NaCl, 15mmol/L NaN3, PH7.2
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09238
Dilution Range	WB:1:500-1:2000, IHC-P:1:20-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse kidney ti...



Immunohistochemical staining of human ki...



Immunofluorescence analysis of HepG2 cel...