


Product Datasheet

60S ribosomal protein L36a-like antibody (HRP) (orb239135)

Description	60S ribosomal protein L36a-like antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human 60S ribosomal protein L36a-like protein (1-106AA)
Target	RPL36AL
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q969Q0
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T whole cell...