


Product Datasheet

60S ribosomal protein L36a-like antibody (FITC) (orb239137)

Description	60S ribosomal protein L36a-like antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human 60S ribosomal protein L36a-like protein (1-106AA)
Target	RPL36AL
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q969Q0
Dilution Range	Not yet tested in other applications.
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1): EC109...