



Product Datasheet

DNA-binding protein inhibitor ID-2 antibody (orb239175)


Description	DNA-binding protein inhibitor ID-2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human DNA-binding protein inhibitor ID-2 protein (1-134AA)
Target	ID2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q02363
Dilution Range	IHC-P:1:1000-1:2000, IF:1:200-1:500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cel...



Immunofluorescence analysis of HepG2 cel...



Immunohistochemical staining of human br...