


Product Datasheet

HCLS1-associated protein X-1 antibody (Biotin) (orb239216)

| | |
|----------------------------|--|
| Description | HCLS1-associated protein X-1 antibody (Biotin) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Biotin |
| Tested Applications | ELISA |
| Immunogen | Recombinant Human HCLS1-associated protein X-1 protein (2-250AA) |
| Target | HAX1 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | 000165 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of (Lane 1): EC109...