



Product Datasheet

Lys-63-specific deubiquitinase BRCC36 antibody (orb239247)


Description	Lys-63-specific deubiquitinase BRCC36 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Lys-63-specific deubiquitinase BRCC36 protein (14-316AA)
Target	BRCC3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P46736
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human re...



Immunohistochemical staining of human ki...



Western blot analysis of HeLa whole cell...