

Product Datasheet

papillomavirus type 18 Regulatory protein E2 antibody (Biotin) (orb239261)

Description	papillomavirus type 18 Regulatory protein E2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human papillomavirus type 18 Regulatory protein E2 protein (1- 365AA)
Target	E2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P06790
Expiration Date	12 months from date of receipt.