

Product Datasheet

papillomavirus type 18 Regulatory protein E2 antibody (HRP) (orb239263)

| | |
|----------------------------|----------------------------------------------------------------------------------|
| Description | papillomavirus type 18 Regulatory protein E2 antibody (HRP) |
| Species/Host | Rabbit |
| Reactivity | Virus |
| Conjugation | HRP |
| Tested Applications | ELISA |
| Immunogen | Recombinant Human papillomavirus type 18 Regulatory protein E2 protein (1-365AA) |
| Target | E2 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P06790 |
| Expiration Date | 12 months from date of receipt. |