

Product Datasheet

Replication protein A 14 kDa subunit antibody (orb239285)

Description Replication protein A 14 kDa subunit antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, WB

Immunogen Recombinant Human Replication protein A 14 kDa subunit protein (1-119AA)

Target RPA3

Preservatives Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.


Note For research use only

Isotype IgG

Clonality Polyclonal

Uniprot ID [P35244](#)

Expiration Date 12 months from date of receipt.



Western blot analysis of LO2(Lane 1), He...