



Product Datasheet

RNA-binding motif protein antibody (orb239304)


| | |
|----------------------------|---|
| Description | RNA-binding motif protein antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF, IHC, WB |
| Immunogen | Recombinant Human RNA-binding motif protein, Y chromosome, family 1 member A1 protein (1-459AA) |
| Target | RBMY1A1 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | PODJD3 |
| Dilution Range | WB:1:500-1:5000, IHC-P:1:20-1:200, IF:1:50-1:200 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of MCF7(Lane 1), 2...



Immunohistochemical staining of human



Western blot analysis of Rat spleen tiss...