


Product Datasheet

C-C motif chemokine 5 antibody (Biotin) (orb239439)

Description	C-C motif chemokine 5 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human C-C motif chemokine 5 protein (24-91AA)
Target	CCL5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P13501
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1): EC109...