

Product Datasheet

C-X-C motif chemokine 10 antibody (orb239461)

| | |
|----------------------------|---|
| Description | C-X-C motif chemokine 10 antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA |
| Immunogen | Recombinant Human C-X-C motif chemokine 10 protein (22-98AA) |
| Target | CXCL10 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P02778 |
| Expiration Date | 12 months from date of receipt. |