

Product Datasheet

High affinity immunoglobulin epsilon R subunit gamma antibody (Biotin) (orb239510)

Description	High affinity immunoglobulin epsilon R subunit gamma antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human High affinity immunoglobulin epsilon receptor subunit gamma protein (19-86AA)
Target	FCER1G
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P30273
Expiration Date	12 months from date of receipt.