



Product Datasheet

Interferon-induced protein with tetratricopeptide repeats 3 antibody (orb239582)


Description	Interferon-induced protein with tetratricopeptide repeats 3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human Interferon-induced protein with tetratricopeptide repeats 3 protein (1-490AA)
Target	IFIT3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14879
Dilution Range	WB:1:1000-1:5000, IHC-P:1:20-1:200, IF:1:50-1:500, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ki...



Immunohistochemical staining of human CO...



Immunoprecipitation analysis of HepG2 wh...