



Product Datasheet

Interferon-induced transprotein 1 antibody (orb239586)

Description	Interferon-induced transprotein 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Interferon-induced transmembrane protein 1 protein (1-36AA)
Target	IFITM1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P13164
Dilution Range	WB:1:1000-1:5000, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of HeLa cells...



Western blot analysis of EC109 whole ce...