



Product Datasheet

NFkB p65 antibody (orb239734)


Description	NFkB p65 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IF, IHC, WB
Immunogen	Recombinant Human Transcription factor p65 protein (1-210AA)
Target	RELA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q04206
Dilution Range	WB:1:500-1:5000, IHC-P:1:500-1:1000, IF:1:200-1:500
Expiration Date	12 months from date of receipt.



Immunoprecipitation analysis of rabbit l...



Immunofluorescence analysis of HepG2 cel...



Immunohistochemical staining of human lu...