

Product Datasheet

Aldehyde oxidase antibody (HRP) (orb239839)

Description	Aldehyde oxidase antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Aldehyde oxidase protein (236-421AA)
Target	Sporamin A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q06278
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.