

Product Datasheet

Fucose-1-phosphate guanylyltransferase antibody (orb239928)

Description	Fucose-1-phosphate guanylyltransferase antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Carcinoembryonic antigen-related cell adhesion molecule 4 protein (36-155AA)
Target	CEACAM4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O75871
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.