



Product Datasheet

Galactose-1-phosphate uridylyltransferase antibody (orb239932)


Description	Galactose-1-phosphate uridylyltransferase antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human Galactose-1-phosphate uridylyltransferase protein (1-379AA)
Target	GALT
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07902
Dilution Range	WB:1:500-1:5000, IHC-P:1:20-1:200, IF:1:50-1:200, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse small int...



Immunoprecipitation analysis of HepG2 wh...



Immunofluorescence analysis of HeLa cell...