

Product Datasheet

O-phosphoseryl-tRNA(Sec) selenium transferase antibody (Biotin) (orb240038)

Description	O-phosphoseryl-tRNA(Sec) selenium transferase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human O-phosphoseryl-tRNA(Sec) selenium transferase protein (5-451AA)
Target	SEPSECS
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HD40
Expiration Date	12 months from date of receipt.