
Product Datasheet

**Activator of 90 kDa heat shock protein ATPase
homolog 1 antibody (orb240199)**

Description

Activator of 90 kDa heat shock

antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ELISA, IHC, WB

Immunogen

Recombinant Human Activator of 90 kDa heat shock protein ATPase homolog 1 protein (1-338AA)

Target

AHS1

PreservativesPreservative: 0.03% Proclin 300.
Constituents: 50% Glycerol**Form/Appearance**

Liquid

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID**O95433****Dilution Range**

WB:1:1000-1:5000, IHC-P:1:20-1:200

Expiration Date

12 months from date of receipt.