


Product Datasheet

Activator of 90 kDa heat shock protein ATPase homolog 1 antibody (Biotin) (orb240200)

Description	Activator of 90 kDa heat shock protein ATPase homolog 1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Activator of 90 kDa heat shock protein ATPase homolog 1 protein (1-338AA)
Target	AHSA1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	095433
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1): EC109...