



Product Datasheet

ADRBK1 antibody (orb240219)


Description	ADRBK1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Beta-adrenergic receptor kinase 1 protein (2-221AA)
Target	GRK2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P25098
Dilution Range	IHC-P:1:500-1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human prostate tissue.



Immunohistochemical staining of human spleen tissue.



Immunohistochemical staining of human tonsil tissue.