



Product Datasheet

Rho guanine nucleotide exchange factor 7 antibody (orb240340)


Description	Rho guanine nucleotide exchange factor 7 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP
Immunogen	Recombinant Human Rho guanine nucleotide exchange factor 7 protein (1-250AA)
Target	ARHGEF7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q14155
Dilution Range	IF:1:200-1:500, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T whole cell...



Immunoprecipitation analysis of HeLa who...



Immunofluorescence analysis of HeLa cell...