


Product Datasheet

Arabidopsis thaliana L-ascorbate peroxidase 2, cytosolic antibody (orb240390)

Description	Arabidopsis thaliana L-ascorbate peroxidase 2, cytosolic antibody
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Arabidopsis thaliana L-ascorbate peroxidase 2, cytosolic protein (4-250AA)
Target	APX2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q1PER6
Expiration Date	12 months from date of receipt.



Western blot analysis of Arabidopsis tha...