

Product Datasheet

Arabidopsis thaliana L-ascorbate peroxidase 2, cytosolic antibody (FITC) (orb240392)

Description	Arabidopsis thaliana L-ascorbate peroxidase 2, cytosolic antibody (FITC)
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	FITC
Immunogen	Recombinant Arabidopsis thaliana L-ascorbate peroxidase 2, cytosolic protein (4-250AA)
Target	APX2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q1PER6
Expiration Date	12 months from date of receipt.