


## Product Datasheet

# Arabidopsis thaliana L-ascorbate peroxidase 2, cytosolic antibody (HRP) (orb240393)

<b>Description</b>	Arabidopsis thaliana L-ascorbate peroxidase 2, cytosolic antibody (HRP)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	A. thaliana
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Recombinant Arabidopsis thaliana L-ascorbate peroxidase 2, cytosolic protein (4-250AA)
<b>Target</b>	APX2
<b>Preservatives</b>	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q1PER6</a>
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of Arabidopsis tha...