

Product Datasheet

Arabidopsis thaliana Protein ABSCISIC ACID-INSENSITIVE 5 antibody (orb240402)

Description	Arabidopsis thaliana Protein ABSCISIC ACID-INSENSITIVE 5 antibody
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Arabidopsis thaliana Protein ABSCISIC ACID-INSENSITIVE 5 protein (1-442AA)
Target	ABI5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9SJN0
Expiration Date	12 months from date of receipt.