

Product Datasheet

Aspergillus niger Probable alpha-galactosidase B antibody (Biotin) (orb240412)

Description	Aspergillus niger Probable alpha-galactosidase B antibody (Biotin)
Species/Host	Rabbit
Reactivity	Fungi
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Aspergillus niger Probable alpha-galactosidase B protein (21-443AA)
Target	agIB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y865
Expiration Date	12 months from date of receipt.