

Product Datasheet

Glyceraldehyde-3-phosphate dehydrogenase antibody (Biotin) (orb240503)

Description	Glyceraldehyde-3-phosphate dehydrogenase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Gallus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Chicken Glyceraldehyde-3-phosphate dehydrogenase protein (2-333AA)
Target	GAPDH
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P00356
Expiration Date	12 months from date of receipt.