

Product Datasheet

Chimpanzee Lymphotoxin-beta antibody (Biotin) (orb240519)

Description	Chimpanzee Lymphotoxin-beta antibody (Biotin)
Species/Host	Rabbit
Reactivity	Primate
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Chimpanzee Lymphotoxin-beta protein (49-244AA)
Target	LTB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q862Z7
Expiration Date	12 months from date of receipt.