

Product Datasheet

FAM127A antibody (FITC) (orb24065)

Description	FAM127A antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Retrotransposon Gag-like protein 8B protein (1-113AA)
Target	RTL8C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A6ZKI3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.